

POST-VASECTOMY SEMEN SAMPLE (PVSA)

Use the **Return Box** for your post-vasectomy semen sample (PVSA)

Send sample **12 weeks** and **20 ejaculations, AFTER** vasectomy. Abstain from ejaculations 2 days before collecting the sample.

Box is labeled be returned to our office. Please make sure you apply **CORRECT POSTAGE**, otherwise it will NOT be delivered. It's better to take it to the Post Office (get a **tracking number**)



Contents of the box:

1. Small **vial** for the sample (it is coated with antibiotic, do NOT wash)
2. Plastic **pipette** (to transfer the semen sample into the vial)
3. Plastic laboratory **bag** with label.

How to SUBMIT your sample?

1. **FILL OUT** the information on the label of the plastic bag.
2. **COLLECT** semen in a clean recipient (ie. clean medicine bottle)
3. **TRANSFER** semen to the vial using the plastic pipette.
4. Close the vial and place it **INSIDE** the plastic bag.
5. Place **bag inside** the Return Box, close it and seal it with tape.
6. Take Return Box **to post office** and **PAY ADEQUATE** postage.

* If you **LOSE** your box, we can send you another one for \$20.